


# Sonic

Aurelia® Sonic® Nitrile Indigo Blue  
Nitrile Examination Gloves



**Indigo Blue Color**  
**300/box Saver Pack**

- New ultra comfort formulation
- 2.2 mil thickness for improved tactile sensitivity
- Finger textured finish
- Cost saving & Space saving: 300 gloves per box
- Soothing indigo blue color



**AURELIA®**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Sonic

Aurelia is proud to introduce our latest innovation in personal protection, Aurelia Sonic® Nitrile Gloves. Featuring our revolutionary ultra-thin formulation Sonic® will adapt so seamlessly to your hands that you'll be working faster than the speed of sound! Each box of Aurelia Sonic® contains 300 machine packed soothing indigo blue colored gloves. We guarantee Sonic® will be the next big BOOM for healthcare professionals! For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Sonic™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	680
Small	31.1	680
Medium	31.0	680
Large	31.2	680
Extra Large	31.3	680
ASTM D6319	14 MPa minimum	500% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	620
Small	30.0	620
Medium	30.0	620
Large	30.0	620
Extra Large	30.0	620
ASTM D6319	14 MPa minimum	620% minimum

### DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Small</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Medium</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Large</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum
<b>Extra Large</b> Aurelia Sonic ASTM D6319	9.4 / 240 8.7 / 220 minimum	2.2 / 0.06 2.0 / 0.05 minimum	3.2 / 0.08 2.0 / 0.05 minimum

### Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 300 pieces per dispenser box, 10 dispenser boxes per shipping carton.

\*Extra Large 270 gloves/ per box

### ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large*
93775	93776	93777	93778	93779

Authorized Aurelia Distributor:



**AURELIA**  
Sign of Comfort™

AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aureliagloves.com](mailto:info@aureliagloves.com) | Facebook: [facebook.com/aureliagloves](https://www.facebook.com/aureliagloves)