


BoldMax

Aurelia® BoldMAX® Black Nitrile
Powder Free Examination Gloves



**Durable Black Nitrile
Patented "MaxGrip" Texture**



Tested and Approved for Fentanyl
and Chemotherapy Drug Resistance

- 6 mil thickness, Patented "MaxGrip" Texture
- Excellent chemical resistance and tensile strength
- Tested and Approved for Fentanyl & Chemo drug resistance per ASTM D6978-05 Standard



AURELIA®
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

BoldMAX

Product Description

BoldMAX™ professional grade black nitrile gloves offers industrial strength and medical grade quality, precision and protection in a stylish matte black finish and are ideal for demanding applications such as automotive, emergency medical and light industrial.



Aurelia® BoldMAX® Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	26.1	550
Small	26.1	550
Medium	26.5	550
Large	26.1	595
Extra Large	26.2	580
XX-Large	26.2	580
ASTM D6319-00	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	24.0	550
Small	24.0	550
Medium	24.0	555
Large	24.0	585
Extra Large	23.7	570
XX-Large	23.7	570
ASTM D6319-00	14 MPa minimum	500% minimum

DIMENSIONES

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum
Small Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum
Large Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum
XX-Large Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	6.7 / 0.13 3.2 / 0.08 minimum	7.9 / 0.15 3.2 / 0.08 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

Medium	Large	Extra Large	XX-Large
97897	97898	97899	97890

Authorized Aurelia Distributor:



AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
 Web: www.aureliagloves.com | Email: info@aureliagloves.com | Facebook: www.facebook.com/aureliagloves